

Accessory Kit for the PerkinElmer DMA 8000

www.lab-ware.co.uk



Ruston Services have developed a small accessory kit of additional items which help to make the PerkinElmer DMA 8000 system even better. The kit includes a stainless steel tray with feet, which provides a waterproof tray to catch melted ice when using the liquid nitrogen cooling system, and also protects your table from liquid nitrogen drips which can damage plastic coatings on laboratory benches. A further benefit is that it can be used to keep the small measuring system parts from falling onto the floor, and is immune to damage if hot parts fall onto it.

Two nut spinners are included for the two sizes of nut which are used with the DMA 8000 measuring systems – 5mm and 5.5mm. These allow easy installation or removal of the nuts on the various threaded legs which comprise the DMA 8000 measuring systems, and help to minimise risk to the DMA drive motor.

A 5mm combination spanner is included with a ring at one end, and open jaw at the other. This is very helpful when installing samples in the tension measuring system.

A USB cooling fan is included, which speeds up cooling of the DMA measuring system and furnace when these are opened at the end of a test. By using the fan, cooling times can be reduced to below 15 minutes from 200°C, giving a faster turnaround between sample tests.

Accessory Kit for the PerkinElmer DMA 8000

www.lab-ware.co.uk

Specifications

Included	Stainless steel tray with rubber feet. Measures 400mm x 200mm x 35mm*
	Strong serrated blade scissors
	Nut spinner for 5mm nuts
	Nut spinner for 5.5mm nuts
	5mm combination spanner with 5mm ring at one end, and 5mm open jaw at the other. Used with the DMA 8000 tension measuring system
	96mm diameter USB Cooling fan with 1m cable to speed up cooling after a test has completed. May be used with USB A extension cables (not included)
	30cm stainless steel safety ruler for cutting samples manually
	Spare stainless steel M2.5 nuts for tension system - pack of 10
	Stainless steel M2.5 x 0.7mm spacers for tension system - pack of 20

* Specification subject to change without notice

Ordering Information

Quotation Requests	Please email your requirements to sales@rustons.org.uk , stating what you would like to order Quotations are normally valid for 60 days.
Ordering	Please send your official order number to the address below or by email to sales@rustons.org.uk If you are outside the UK and your organisation is VAT registered, please ensure that your VAT number is clearly visible on your order.
Delivery	Delivery is normally 2 - 4 weeks from receipt of order. We may be able to provide faster delivery in some circumstances